

Analysis of Intact Proteins on a Thermo Scientific Accucore 150-C4 150 Å Pore  
Diameter NanoLC Column

Application #379

24 August 2012

Copyright © 2010 Thermo Fisher Scientific Inc.

Exported by Thermo Fisher Scientific Applications Library

# Analysis of Intact Proteins on a Thermo Scientific Accucore 150-C4 150 Å Pore Diameter NanoLC Column

## General Information

Market: Biopharmaceutical

Matrix: Standard

Instrument type: UHPLC

## Description

Accucore HPLC columns use Core Enhanced Technology to facilitate fast and high efficiency separations. The 2.6 µm diameter particles have a solid core and a porous outer layer. The optimized phase bonding creates a series of high coverage, robust phases. The tightly controlled 2.6µm diameter of accucore particles results in much lower backpressures than typically seen with 2 µm materials. For the analysis of large biomolecules the Accucore pore size has been further optimized and a C4 phase with reduced hydrophobic retention has been prepared. This 150 Å pore size enables the effective analysis of molecules unable to penetrate into smaller diameter pores and the low hydrophobicity C4 phase results in protein separation by hydrophobicity.

## Method Details

### Instrument parameters

Instrument Parameter	Value
Run Time Length	20.003 min
Wavelength 1	214.0 nm
Wavelength 2	281.0 nm
Mobile_Phase_A	0.1 % formic acid in water
Mobile_Phase_B	0.1 % formic acid in acetonitrile
Flow_Rate	300 nL/min
Backpressure_(100%_aqueous)	204 bar
Run_time	15 minutes + equilibration time
Column_Temperature	40 °C
Injection_details	0.25 µL 2 pmol/µL solution
UV_detector_wavelength	214 nm

## Gradient Details

Ret. Time [min]	Flow [ml/min]	%B [%]	%C [%]	%D [%]
0.000	0.000	0.0	0.0	0.0
1.000	0.000	30.0	0.0	0.0
11.000	0.000	60.0	0.0	0.0
12.000	0.000	95.0	0.0	0.0
15.000	0.000	95.0	0.0	0.0

## Column Details


column_A	
Manufacturer	Thermo Fisher Scientific
Model	Accucore 150-C4
Diameter	75
Length	150
Particle Size	2.6
Packing Material	150-C4

## System information

Instrument Type	UHPLC
UHPLC	Generic
UV Detector model	VWD-3400


## Results

### Channel UV\_VIS\_1



No	Peak_Name	*Compound _Class	Retention_Ti me	Peak_Area	Peak_Area_p c	Peak_Height	Peak_Height _pc	Plates_(USP)	Resolution_ (USP)	Tailing_Facto r_(USP)
1	cytochrome C	Protein	6.940	1.143	9.95	10.162	16.29	24095	3.24	1.18
2	Insulin	Protein	7.403	1.008	8.78	12.712	20.38	76490	7.10	0.99
3	Myoglobin	Protein	8.823	1.720	14.98	8.855	14.19	14490	2.41	2.18
4	Carbonic Anhydrase	Protein	9.590	3.804	33.13	17.069	27.36	12507	3.86	2.24
5	Ovalbumin	Protein	11.023	3.808	33.16	13.590	21.78	12176	n.a.	1.81

Channel UV\_VIS\_2



No	Peak_Name	*Compound _Class	Retention_Ti me	Peak_Area	Peak_Area_p c	Peak_Height	Peak_Height _pc	Plates_(USP)	Resolution_( USP)	Tailing_Facto r_(USP)
1	Carbonic Anhydrase	Protein	9.590	0.335	60.81	1.590	62.38	23079	4.34	2.03
2	Ovalbumin	Protein	11.053	0.216	39.19	0.959	37.62	11031	n.a.	1.12

## Appendix

The application can be accessed at <http://dlibrary.dionex.com/Public/View.aspx?ApplicationID=379>

### Available Downloads

Filename	Size(bytes)
Analysis of Intact Proteins.pdf	616672

### Related Information

No related information available.