



#### Antiepileptics

Column: YMC-Pack Pro C18 (5  $\mu$ m, 12 nm, 150 x 4.6 mm)

Eluent: 50 mM NH<sub>4</sub>H<sub>2</sub>PO<sub>4</sub>-(NH<sub>4</sub>)<sub>2</sub>HPO<sub>4</sub> (pH 6.7) / methanol (50/50)

Flow: 0.7 ml/min

Detection: UV at 215 nm, 0.32 AUFS

Temperature: 30 °C

Injection: 8  $\mu$ l (0.035  $\square$  0.7 mg/ml)