

Neue Test- Separation of hydrophobic, polar and basic mixture

Column:	Listed on chromatogram
Dimension:	50 x 2.1 mm ID
Mobile Phase:	MeOH:20 mM phosphate buffer (pH 7.0) = 80:20
Flow Rate:	0.5 ml/min
Temperature:	Ambient
Detection:	UV 254 nm
Sample:	<ol style="list-style-type: none">1. Uracil2. Butylparaben3. Dipropyl phthalate4. Naphthalene5. Acenaphthene6. Amitriptyline

*Eclipse is registered trademarks of Agilent Technologies. Kinetex is registered trademarks of Phenomenex. Acquity is registered trademarks of Waters Corporation. Dikma Technologies is not affiliated with the above companies.

Neue Test- Separation of hydrophobic, polar and basic mixture

