

Antibacterials

Column: 50 x 2.1 mm ID

Mobile phase: MeCN + 0.1% HCOOH:H₂O +0.1% HCOOH = 20:80

Flow rate: 0.3 mL/min

Detection: UV 254 nm

Sample: 1. Sulfanilamide

2. Carbadox

3. Sulfamerazine

4. Sulfamethazine

5. Sulfamethoxypyridazine

6. Furazolidone

7. Sulfamethoxazole

8. Sulfisoxazole

9. Orolinic acid

10. Sulfadimethoxine

