

Cold Medicine

Columns:	Listed on chromatograms
Dimensions:	150×4.6mm ID, 5µm(all columns)
Mobile Phase:	ACN+0.1%TFA: H₂O+0.1%TFA=35:65
Flow Rate:	1.0ml/min
Temperature:	Ambient
Detection:	UV 220nm
Sample:	1. Doxylamine
	2. Acetanilide
	3. Aspirin
	4. Dextromethorphan
	5. Diphenhydramine

*ZORBAX,HC,TC,Eclipse,Hypersil are registered trademarks of Agilent Technologies.Gemini,Luna are registered trademarks of Phenomenex.Symmetry,XBridge are registered trademarks of Waters Corporation.Inertsil is registered trademark of GL Sciences Inc.Kromasil is registered trademark of Eka Chemicals.BECKMAN Coulter is registered trademark of Beckman Coulter, Inc. SMT is registered trademark of Analtech Inc.Dikma Technologies is not affiliated with the above companies.









