

# Multiple Injections

<b>Columns:</b>	<b>Spursil C18</b>
<b>Dimensions:</b>	<b>150 × 4.6mm 5 μm</b>
<b>Mobile Phase:</b>	<b>MeOH: 0.1% TFA in H<sub>2</sub>O =30:70</b>
<b>Flow Rate:</b>	<b>1.0ml/min</b>
<b>Temperature:</b>	<b>Ambient</b>
<b>Detection:</b>	<b>UV 230nm</b>
<b>Sample:</b>	<b>1. Ceftazidime</b>
	<b>2. Cefadroxil</b>
	<b>3. Cefuroxime</b>
	<b>4. Cefazoline</b>
	<b>5. Cefaclor</b>
	<b>6. Cephalexin</b>

