



test parameters	column: Astec (R,R) P-CAP or (S,S) P-CAP, 25 cm x 4.6 mm I.D., 5 µm particles (31024AST, 33024AST)
	mobile phase: (A) acetonitrile; (B) methanol; (97:3, A:B)
	flow rate: 1 mL/min
	column temp.: 25 °C
	detector: UV, 254 nm
	injection: 10 µL
	sample: (R,R)- and (S,S)-DACH-ACR (N-(2- acryloylamino-(1R,2R)-cyclohexyl)- acrylamide)