


1. Thyroglobulin, 669,000 Da
2. Alcohol dehydrogenase, 150,000 Da
3. Albumin, 66,000 Da
4. Carbonic anhydrase, 29,000 Da
5. Aprotinin, 6,500 Da
6. Cyanocobalamin, 1,355 Da
7. p-Aminobenzoic acid, 137 Da



test parameters	column: Discovery BIO GFC 150, 30 cm x 7.8 mm I.D., 5 µm particles (567300-U)
	mobile phase: 150 mM potassium phosphate, pH 7.0
	flow rate: 1.0 mL/min
	column temp.: 25 °C
	detector: UV, 280 nm
	injection: 10 µL
	sample: 1-10 g/L in buffer