



Sample : 20 μ L

1. 5'CMP 20 μ g/mL
2. 5'UMP 20 μ g/mL
3. 5'TMP 20 μ g/mL
4. 3'CMP 20 μ g/mL
5. 5'ADP 100 μ g/mL
6. 3'UMP 20 μ g/mL
7. 5'GMP 20 μ g/mL
8. 2'AMP 20 μ g/mL
9. 3'AMP 20 μ g/mL
10. 3'GMP 20 μ g/mL

Column	: Shodex AXpak WA-624 (6.0mmID*150mm)
Eluent	: 0.35M CH ₃ COOH aq./0.35M CH ₃ COONH ₄ aq.=240/100
Flow rate	: 1.0mL/min
Detector	: UV(260nm)
Column temp.	: 60°C