
Catecholamine Analysis (II)

ACE®
Ultra-inert
UHPLC & HPLC Columns

Application #AN3910

Conditions

Column:	ACE 5 AQ
Dimensions:	150 x 4.6 mm
Part Number:	ACE-126-1546
Mobile Phase:	50 mM KH_2PO_4 pH 3.0 in H_2O
Flow Rate:	1 mL/min
Injection:	2 µL
Temperature:	Ambient
Detection:	UV, 210 nm



1. Noradrenaline
2. Adrenaline
3. L-DOPA
4. Dopamine
5. L-Tyrosine
6. VMA (Vanillylmandelic acid)

© 2015 ACT Ltd. All rights reserved. All products are available worldwide.