Extract of Ginseng Stems and Leaves

ACE®
Ultra-inert
UHPLC & HPLC Columns


Application #AN6700

Conditions

Columns:	ACE UltraCore 2.5 SuperC18
Dimensions:	50 x 4.6 mm
Part Numbers:	CORE-25A-0546U
Mobile Phase:	A: H_2O B: MeCN

Time (mins)	%B
0	19
6	19
10	29
13	29
18	40
25	19

Flow Rate:	1.3 mL/min
Injection:	10 µL
Temperature:	35 °C
Detection:	UV, 203 nm



1. Ginsenoside Rg1
2. Ginsenoside Re
3. Ginsenoside Rf
4. Ginsenoside Rb1
5. Ginsenoside Rc
6. Ginsenoside Rb2
7. Ginsenoside Rd

© 2018 ACT Ltd. All rights reserved. All products are available worldwide.