
Determination of Paeonol in Mingmu Dihuang Pills

ACE®
Ultra-inert
UHPLC & HPLC Columns


Application #AN6550

Conditions

Columns:	ACE Excel 5 C18
Dimensions:	250 x 4.6 mm
Part Numbers:	EXL-121-2546U
Mobile Phase:	MeOH/H_2O (48:52 v/v)
Flow Rate:	1.0 mL/min
Injection:	10 µL
Temperature:	Ambient
Detection:	UV, 274 nm
Method:	ChP.2015, Vol. 1



1. Paeonol



© 2018 ACT Ltd. All rights reserved. All products are available worldwide.