
Determination of Isoimperatorin and Osthol in Tianma pills

ACE®
Ultra-inert
UHPLC & HPLC Columns

Application #AN6030

Conditions

Column:	ACE Excel 5 C18-PFP
Dimensions:	250 x 4.6 mm
Part Number:	EXL-1210-2546U
Mobile Phase:	MeCN/0.1 % phosphoric acid (aq.) (58:42 v/v)
Flow Rate:	1.0 mL/min
Injection:	10 µL
Temperature:	45 °C
Detection:	UV, 320 nm
Monograph Method:	ChP.2015, Vol. 1



1. Osthol

2. Isoimperatorin

© 2018 ACT Ltd. All rights reserved. All products are available worldwide.