
Aspirin and Related Substances (III)


ACE®
Ultra-inert
UHPLC & HPLC Columns

Application #AN4740


Conditions

Column:	ACE Excel 1.7 SuperC18
Dimensions:	50 x 3.0 mm
Part Number:	EXL-1711-0503U
Mobile Phase:	MeCN/MeOH/H_2O/H_3PO_4 (35:5:60:0.2 v/v/v/v)
Flow Rate:	1.25 mL/min
Injection:	0.7 µL
Temperature:	40 °C
Detection:	UV, 237 nm
Instrument:	Chromaster UltraRs







1. 4-Hydroxybenzoic acid




2. 4-Hydroxyisophthalic acid




3. Acetylsalicylic acid (Aspirin)



4. Salicylic acid



5. Acetylsalicylsalicylic acid



6. Salsalate

© 2016 ACT Ltd. All rights reserved. All products are available worldwide.