
Chromatography Solutions

Application note #7090

Separation of Psoralen and Isopsoralen




Avantor® ACE®


Method Details

CONDITIONS

Column:	Avantor® ACE® Excel 5 C18-PFP
Particle Size:	5 µm
Dimensions:	250 x 4.6 mm
Mobile Phase:	H_2O/MeOH (40:60 v/v)
Flow Rate:	1.0 mL/min
Injection:	10 µL
Detection:	UV, 240 nm



Psoralen



Isopsoralen

ORDERING TABLE

Product	Details	Size	Part Number
Avantor® ACE® Excel 5 C18-PFP	HPLC Column	250 x 4.6 mm	EXL-1210-2546U