
Antihistamines and Expectorants – Mobile Phase Effects

ACE®
Ultra-inert
UHPLC & HPLC Columns


Application #AN3960

Conditions

Column:	ACE 5 C18
Dimensions:	150 x 4.6 mm
Part Number:	ACE-121-1546
Mobile Phase:	(a) MeOH/50 mM KH_2PO_4 pH 3.0 in H_2O (45:55 v/v)
	(b) MeCN/50 mM KH_2PO_4 pH 3.0 in H_2O (28:72 v/v)
Flow Rate:	1 mL/min
Temperature:	Ambient
Detection:	UV, 220 nm



1. Maleic acid
2. Salicylamide
3. Guaifenesin
4. Guaiacol
5. Chlorpheniramine maleate
6. Dextromethorphan



© 2015 ACT Ltd. All rights reserved. All products are available worldwide.