
Sugars – Cola Vs Diet Cola

ACE®
Ultra-inert
UHPLC & HPLC Columns

Application #AN4150

Conditions

Column:	ACE Excel 3 NH_2
Dimensions:	150 x 4.6 mm
Part Number:	EXL-1114-1546U
Mobile Phase:	MeCN/H_2O (75:25 v/v)
Flow Rate:	1 mL/min
Injection:	10 µL
Temperature:	35 °C
Detection:	RI, 35 °C
Sample:	100 µL drink in 900 µL mobile phase
System:	Chromaster 600



1. Fructose
2. Glucose
3. Sucrose

© 2015 ACT Ltd. All rights reserved. All products are available worldwide.