
# USP Method for Analysis of Povidone Impurity B (2-Pyrrolidone)

ACE®
Ultra-inert
UHPLC & HPLC Columns

**Application #AN6840**

**Conditions**

| | |
|---|---|
| Columns: | ACE 5 C18 |
| Dimensions: | 150 x 4.6 mm |
| Part Number: | ACE-121-1546 |
| Mobile Phase: | 5:95 MeOH/$H_2O$ |
| Flow Rate: | 0.8 mL/min |
| Injection: | 5 µL |
| Temperature: | 40 °C |
| Detection: | UV, 205 nm |







2-Pyrrolidone (Impurity B)

Always check the latest and official method information from the relevant pharmacopoeia prior to analysis.

© 2018 ACT Ltd. All rights reserved. All products are available worldwide.