
Analysis of Povidone Impurities


ACE®
Ultra-inert
UHPLC & HPLC Columns

Application #AN6820


Conditions

Columns:	ACE 5 C18
Dimensions:	150 x 4.6 mm
Part Number:	ACE-121-1546
Mobile Phase:	10:90 MeCN/H_2O
Flow Rate:	1.0 mL/min
Injection:	5 µL
Temperature:	40 °C
Detection:	UV, 205 nm





1. 2-pyrrolidone (Impurity B)



2. Vinylpyrrolidinone (Impurity A)

© 2018 ACT Ltd. All rights reserved. All products are available worldwide.