
Chromatography Solutions

Application note #1030

Corticosteroids by LC-MS/MS



Avantor® ACE®

Method Details

CONDITIONS

Column: ACE 3 C18-PFP

Dimensions: 150 x 2.1 mm

Part Number: ACE-1110-1502

Mobile Phase: A: 0.1% formic acid in H_2O

B: MeCN/0.1% formic acid in H_2O (90:10 v/v)

Gradient:

Time (mins)	%B
0.0	30
14.0	50
17.0	95
20.0	30

Flow Rate: 0.3 mL/min

Injection: 25 µL

Temperature: 15 °C

Detection: Turbospray, MRM

ORDERING TABLE

Product	Details	Size	Part number
Avantor® ACE® 3 C18-PFP	HPLC Column	150 x 2.1 mm	ACE-1110-1502

Reproduced with permission of Laboratorio Salud Publica de Burgos, Spain

vwr.com/ace
chromsupport@avantorsciences.com