Pharmaceutically Relevant Compounds (I)

ACE®
Ultra-inert
UHPLC & HPLC Columns

Application #AN1580

Conditions


Column:	ACE 5 C18-AR
Dimensions:	150 x 4.6 mm
Part Number:	ACE-129-1546
Mobile Phase:	A: 5 mM formic acid in H_2O B: 5 mM formic acid in MeOH

Time (mins)	%B
0.0	3
5.0	100
5.5	100
6.0	3

Post time 2.5 minutes

Flow Rate:	1.5 mL/min
Injection:	5 µL
Temperature:	40 °C
Detection:	UV, 254 nm

1. Paracetamol
2. Hydrochlorothiazide
3. Methylphenylsulfoxide
4. Methylphenylsulfone
5. Aspirin
6. Phenacetin
7. 1,3-Dinitrobenzene
8. 1,2,4-Trimethoxybenzene
9. Ethylbenzoate
10. Nimesulide
11. Ibuprofen
12. Indomethacin
13. Mefenamic acid



© 2015 ACT Ltd. All rights reserved. All products are available worldwide.