
Determination of Loganin in Mingmu Dihuang Pills

ACE®
Ultra-inert
UHPLC & HPLC Columns


Application #AN6540

Conditions

Columns:	ACE Excel 5 C18
Dimensions:	250 x 4.6 mm
Part Numbers:	EXL-121-2546U
Mobile Phase:	[THF/MeCN/MeOH (1:8:4 v/v/v)]/0.05 % H_3PO_4 solution (8:92 v/v)
Flow Rate:	1.0 mL/min
Injection:	10 µL
Temperature:	Ambient
Detection:	UV, 236 nm
Method:	ChP.2015, Vol. 1



1. Loganin



© 2018 ACT Ltd. All rights reserved. All products are available worldwide.