Pharmaceutically Relevant Compounds and Comparison (II)

ACE® Ultra-inert UHPLC & HPLC Columns


Application #**AN1630**


Conditions

Column:	ACE Excel 3 CN-ES ACE Excel 3 C18-Amide ACE Excel 3 C18-PFP
Dimensions:	100 x 3.0 mm
Part Numbers:	EXL-1113-1003U, EXL-1112-1003U, EXL-1110-1003U
Mobile Phase:	A: 20 mM ammonium formate in H_2O B: 20 mM ammonium formate in MeOH


Time (mins)	%B
0.0	3
5.0	100
6.0	100
6.5	3

Flow Rate:	0.6 mL/min
Temperature:	40 °C
Detection:	UV







1. Hydrochlorothiazide




2. Methylphenylsulfoxide




3. 1,3,5-Trinitrobenzene




4. Myricetin



5. p-Cresol



6. Sulindac



7. Toluene

© 2015 ACT Ltd. All rights reserved. All products are available worldwide.