Ginkgo Biloba – Ultra Resolution

ACE®
Ultra-inert
UHPLC & HPLC Columns


Application #AN2270

Conditions

Column:	ACE Excel 1.7 C18-PFP
Dimensions:	2 x 100 x 3.0 mm (coupled)
Part Number:	2 x EXL-1710-1003U
Mobile Phase:	A: 0.1% formic acid in H_2O
	B: 0.1% formic acid in MeCN

Time (mins)	%B
0.00	5
0.72	5
15.72	50
18.72	100
20.72	100
22.72	5

Flow Rate:	0.8 mL/min
Injection:	2 µL
Temperature:	80 °C
Detection:	UV, 254 nm



1. Ginkgo Biloba

© 2015 ACT Ltd. All rights reserved. All products are available worldwide.