
Determination of Amoxicillin Content

ACE®
Ultra-inert
UHPLC & HPLC Columns


Application #AN6040

Conditions

Column:	ACE Excel 5 SuperC18
Dimensions:	250 x 4.6 mm
Part Number:	EXL-1211-2546U
Mobile Phase:	50 mM KH_2PO_4 (pH 5)/MeCN (97.5:2.5)
Flow Rate:	1.0 mL/min
Injection:	20 µL
Temperature:	Ambient
Detection:	UV, 254 nm
Monograph Method:	ChP. 2015, Vol. 2



1. Amoxicillin



© 2018 ACT Ltd. All rights reserved. All products are available worldwide.