
Analysis of Isoniazid and Related Impurities

ACE® Ultra-inert UHPLC & HPLC Columns


Application #**AN6600**

Conditions

Columns:	ACE Excel 5 C18-AR
Dimensions:	250 x 4.6 mm
Part Numbers:	EXL-129-2546U
Mobile Phase:	MeOH/20 mM Na_2HPO_4 solution (pH 6) (15:85 v/v)
Flow Rate:	1.0 mL/min
Injection:	10 µL
Temperature:	Ambient
Detection:	UV, 262 nm
Method:	ChP.2015, Vol. 2



1. Impurity A
2. Isoniazid
3. Impurity B
4. Impurity C
5. Impurity D
6. Impurity E



© 2018 ACT Ltd. All rights reserved. All products are available worldwide.