Antihistamines

ACE®
Ultra-inert
UHPLC & HPLC Columns


Application #AN1400

Conditions

Column:	ACE UltraCore 2.5 SuperC18
Dimensions:	100 x 3.0 mm
Part Number:	CORE-25A-1003U
Mobile Phase:	A: 20 mM ammonium formate pH 3.0 in H_2O
	B: 20 mM ammonium formate pH 3.0 in MeOH/H_2O (9:1 v/v)

Time (mins)	%B
0.0	0
7.5	100
8.5	100
9.0	0

Flow Rate:	1.2 mL/min
Injection:	2 µL
Temperature:	25 °C
Detection:	UV, 254 nm



1. Pseudoephedrine
2. Scopolamine
3. Doxylamine
4. Chlorpheniramine
5. Triprolidine
6. Diphenhydramine
7. Acrivastine
8. Promethazine
9. Fexofenadine
10. Cetirizine
11. Loratadine

© 2015 ACT Ltd. All rights reserved. All products are available worldwide.