
USP Method for Analysis of Acarbose

ACE®
Ultra-inert
UHPLC & HPLC Columns


Application #AN6380

Conditions

Column:	ACE Excel 5 NH2
Dimensions:	250 x 4.6 mm
Part Number:	EXL-1214-2546U
Mobile phase:	75:25 MeCN/buffer (0.6 g KH_2PO_4 + 0.35 g Na_2HPO_4 in 1000 mL water)
Flow Rate:	2.0 mL/min
Injection:	10 µL
Temperature:	35 °C
Detection:	UV, 210 nm
Monograph method:	USP 37



1. Acarbose



© 2018 ACT Ltd. All rights reserved. All products are available worldwide.