Additives and Intense Sweeteners

ACE®
Ultra-inert
UHPLC & HPLC Columns


Application #AN2950

Conditions

Column:	ACE 5 C18
Dimensions:	250 x 4.0 mm
Part Number:	ACE-121-2504
Mobile Phase:	A: H_2O
	B: MeCN
	C: 1% TFA in H_2O

Time (mins)	%A	%B	%C
0	88	2	10
25	50	40	10
30	30	60	10
35	88	2	10

Flow Rate:	1 mL/min
Temperature:	30 °C
Detection:	ELSD



1. Acesulfame K
2. Theobromine
3. Theophylline
4. Cyclamate
5. Saccharin
6. Caffeine
7. Sucralose
8. Quinine sulphate
9. Aspartame
10. Neohesperidin dihydrochalcone

Reproduced with permission of Durham County Council, UK
© 2015 ACT Ltd. All rights reserved. All products are available worldwide.