
Amoxicillin USP: Assay

ACE®
Ultra-inert
UHPLC & HPLC Columns

Application #AN4800

Conditions

Column:	ACE 10 C18
Dimensions:	250 x 4.6 mm
Part Number:	ACE-131-2546
Mobile Phase:	50 mM monobasic potassium phosphate pH 5.0: MeCN (24:1 v:v)
Flow Rate:	1.5 mL/min
Injection:	10 µL
Temperature:	22 °C
Detection:	UV, 230 nm
Sample:	Amoxicillin Oral Suspension: 1 mg/mL amoxicillin in 50 mM monobasic potassium phosphate pH 5.0
System:	Chromaster 600



Amoxicillin

© 2017 ACT Ltd. All rights reserved. All products are available worldwide.