
Procaine and p-Aminobenzoic Acid Separation


ACE®
Ultra-inert
UHPLC & HPLC Columns

Application #AN1660


Conditions

Column:	ACE 3 C18-PFP
Dimensions:	100 x 4.6 mm
Part Number:	ACE-1110-1046
Mobile Phase:	0.6% acetic acid in H_2O/MeOH (81:19 v/v), adjusted to pH 4.7 with 20% NaOH
Flow Rate:	1 mL/min
Detection:	UV, 279 nm S/N limit: 10
Sample:	Procaine 0.0002 mg/mL, p-Aminobenzoic acid 0.00005 mg/mL







1. p-Aminobenzoic acid



2. Procaine





Reproduced with permission of Huddersfield Royal Infirmary, UK

© 2015 ACT Ltd. All rights reserved. All products are available worldwide.