
Determination of Impurities in Oxaliplatin and Related Substances

ACE®
Ultra-inert
UHPLC & HPLC Columns


Application #AN6470

Conditions

Column:	ACE Excel 5 AQ
Dimensions:	250 x 4.6 mm
Part Number:	EXL-126-2546U
Mobile Phase:	Potassium phosphate buffer (pH 3.0):MeCN (80:20 v/v) (Dissolve 1.36 g KH_2PO_4 + 1 g sodium 1-heptanesulfonate in 1000 mL H_2O, adjust to pH 3.0)
Flow Rate:	2.0 mL/min
Injection:	20 µL
Temperature:	40 °C
Detection:	UV, 215 nm
Monograph:	EP 9.0



1. Impurity B



2. Impurity E

© 2018 ACT Ltd. All rights reserved. All products are available worldwide.