
Chromatography Solutions

Application note #7530

Hydrocortisone Cream USP Assay



1. Hydrocortisone


Avantor® ACE®

Method Details

CONDITIONS

Columns:	Avantor® ACE® 5 C18
Particle Size:	5 µm
Dimensions:	300 x 3.9 mm
Mobile Phase:	H_2O/MeCN (70:30)
Flow Rate:	1.0 mL/min
Injection Volume:	20 µL
Temperature:	25 °C
Detection:	UV, 254 nm
Method:	USP 37

Always check the latest and official method information from the relevant pharmacopoeia prior to analysis.



Hydrocortisone

ORDERING TABLE

Product	Details	Size	Part number
Avantor® ACE® 5 C18	HPLC Column	300 x 3.9 mm	ACE-121-3039

Reproduced with permission of Mumbai Avantor Application lab.

vwr.com/ace
chromsupport@avantorsciences.com