
# Dexamethasone and Related Substances

ACE®
Ultra-inert
UHPLC & HPLC Columns


**Application #AN6370**

**Conditions**


| | |
|---|---|
| Column: | ACE Excel 5 C18 |
| Dimensions: | 250 x 4.6 mm |
| Part Number: | EXL-121-2546U |
| Mobile phase: | 40 mM $KH_2PO_4$ and 24 mM hexylamine in MeCN:$H_2O$ (71:29 v/v) |
| Flow Rate: | 1.0 mL/min |
| Injection: | 20 µL |
| Temperature: | Ambient |
| Detection: | UV, 254 nm |
| Monograph method: | EP 6.0 |




1. Betamethasone




2. Dexamethasone




© 2018 ACT Ltd. All rights reserved. All products are available worldwide.