
Polyamines

ACE®
Ultra-inert
UHPLC & HPLC Columns

Application #AN3740

Conditions

Column:	ACE 5 C18
Dimensions:	150 x 4.6 mm
Part Number:	ACE-121-1546
Mobile Phase:	TRIS buffer pH 7.0/MeOH (10:90 v/v)
Flow Rate:	1.2 mL/min
Detection:	Fluorescence – λ_{ex} 340 nm, λ_{em} 450 nm



1. Putrescine
2. Cadaverine

(as OPA derivatives)

Reproduced with permission of Centre for Chemical Sciences, Royal Holloway University of London, UK
© 2015 ACT Ltd. All rights reserved. All products are available worldwide.