
Sugars Separation

ACE®
Ultra-inert
UHPLC & HPLC Columns

Application #AN4030

Conditions

Column:	ACE Excel 5 NH_2
Dimensions:	150 x 4.6 mm
Part Number:	EXL-1214-1546U
Mobile Phase:	MeCN/H_2O (70:30 v/v)
Flow Rate:	1 mL/min
Injection:	10 µL
Temperature:	35 °C
Detection:	RI, 35 °C
Sample:	5 mg/mL in mobile phase
System:	Chromaster 600



1. Fructose
2. Sucrose
3. Lactose

© 2015 ACT Ltd. All rights reserved. All products are available worldwide.