Biomarker Profiling


ACE®
Ultra-inert
UHPLC & HPLC Columns

Application #AN1990

Conditions

Column:	ACE UltraCore 2.5 SuperC18
Dimensions:	50 x 2.1 mm
Part Number:	CORE-25A-0502U
Mobile Phase:	A: 7% methanol, 2.1 mM tetrabutylammonium bisulphate + 84 mM KH_2PO_4, pH 6 with KOH B: 7% methanol, 2.1 mM tetrabutylammonium bisulphate + 8.4 mM KH_2PO_4, pH 6 with KOH A/B (90:10 v/v)
Flow Rate:	0.4 mL/min
Injection:	2 µL
Temperature:	Ambient
Detection:	UV, 260 nm

1. Deoxyuridine
2. Deoxyguanosine
3. Deoxythymidine
4. Adenosine
5. Deoxyadenosine
6. 2-Fluoro-2'-fluoroadenine arabinoside
7. 2-Fluorodeoxyadenosine
8. 2-Fluoroadenosine
9. 2-Chlorodeoxyadenosine
10. Clofarabine
11. Methylthioadenosine



Reproduced with permission of Department of Medical Biochemistry and Biophysics, Umeå University, Sweden

© 2015 ACT Ltd. All rights reserved. All products are available worldwide.