BSA Tryptic Digest Analysis using Column Coupling

ACE®
Ultra-inert
UHPLC & HPLC Columns


Application #AN6910

Conditions

Column:	ACE UltraCore 2.5 SuperC18
Dimensions:	3 x 150 x 2.1 mm
Part Number:	CORE-25A-1502U
Mobile Phase:	A: 0.05 % TFA in H_2O
	B: 0.05 % TFA in MeCN

Time (mins)	%B
0	10
90	40
120	65
125	95
130	95
132	10
180	10

Flow Rate:	0.21 mL/min
Injection Volume:	20 µL
Temperature:	60 °C
Detection:	214 nm



© 2019 ACT Ltd. All rights reserved. All products are available worldwide.