Pharmaceutically Relevant Molecules (IV)

ACE®
Ultra-inert
UHPLC & HPLC Columns


Application #AN2460

Conditions

Column:	ACE Excel 3 CN-ES
Dimensions:	100 x 2.1 mm
Part Number:	EXL-1113-1002U
Mobile Phase:	A: 0.1% formic acid in H_2O
	B: 0.1% formic acid in MeCN

Time (mins)	%B
0	3
10	100

Flow Rate:	0.6 mL/min
Temperature:	40 °C
Detection:	UV, 210 nm



1. 1,3-Dihydroxybenzene
2. Catechol
3. Hydrochlorothiazide
4. Oxprenolol
5. Salicylic acid
6. Myricetin
7. Piroxicam
8. 1,2-Dinitrobenzene
9. Tolmetin
10. 1-Naphthol
11. Piperine
12. Diflunisal
13. Propylbenzene
14. 1,2,3-Trichlorobenzene

© 2015 ACT Ltd. All rights reserved. All products are available worldwide.