Green Tea Metabolite Profiling by LC-MS

ACE®
Ultra-inert
UHPLC & HPLC Columns


Application #AN2580

Conditions

Column:	ACE Excel 1.7 C18-Amide
Dimensions:	100 x 2.1 mm
Part Number:	EXL-1712-1002U
Mobile Phase:	A: 0.01% formic acid in H_2O B: 0.01% formic acid in MeCN

Time (mins)	%B
0.0	3
2.5	10
8.0	100
8.5	3
10.0	3

Flow Rate:	0.5 mL/min
Detection:	Exactive accurate mass MS system ESI in negative ion mode Analytes between *m/z* 70-800 monitored
Sample:	Metabolites from green tea extracted into cold water by vortexing for 20 mins. Samples filtered prior to injection onto column and modular Accela LC system



Reproduced with permission of School of Pharmacy, University of Nottingham, UK
© 2015 ACT Ltd. All rights reserved. All products are available worldwide.