
Gamma Hydroxybutyric Acid (GHB) and Gamma Butyrolactone (GBL) Separation


ACE®
Ultra-inert
UHPLC & HPLC Columns

Application #AN1500


Conditions

Column:	ACE Excel 2 C18-AR
Dimensions:	100 x 3.0 mm
Part Number:	EXL-109-1003U
Mobile Phase:	20 mM KH_2PO_4 pH 2.5 H_2O/MeCN (98:2 v/v)
Flow Rate:	0.43 mL/min
Injection:	2 µL
Temperature:	30 °C
Detection:	UV, 215 nm





1. GHB



2. GBL

© 2015 ACT Ltd. All rights reserved. All products are available worldwide.