
USP Method for Analysis of Povidone Impurity A (Vinylpyrrolidinone)

ACE®
Ultra-inert
UHPLC & HPLC Columns

Application #**AN6830**

Conditions

Columns:	ACE 5 C18
Dimensions:	150 x 4.6 mm
Part Number:	ACE-121-1546
Mobile Phase:	10:90 MeCN/H_2O
Flow Rate:	1.0 mL/min
Injection:	5 µL
Temperature:	40 °C
Detection:	UV, 235 nm



Vinylpyrrolidinone (Impurity A)

Always check the latest and official method information from the relevant pharmacopoeia prior to analysis.

© 2018 ACT Ltd. All rights reserved. All products are available worldwide.