Green Tea Extract

ACE®
Ultra-inert
UHPLC & HPLC Columns


Application #AN4280

Conditions

Column:	ACE Excel 1.7 SuperC18
Dimensions:	100 x 3.0 mm
Part Number:	EXL-1711-1003U
Mobile Phase:	A: 0.1% formic acid in H_2O B: 0.1% formic acid in MeCN

Time (mins)	%B
0	5
15	20
17	95
18	95
20	5
27	5

Flow Rate:	0.8 mL/min
Injection:	2 µL
Temperature:	80 °C
Detection:	UV, 260 nm
Sample:	Tablet ground to fine powder and extracted with MeCN/H_2O (1:1 v/v) with ultrasonication. Supernatant diluted with H_2O and filtered using Whatman Mini-Uniprep syringeless filter
System:	Chromaster Ultra Rs



© 2015 ACT Ltd. All rights reserved. All products are available worldwide.