Sulfonamides

ACE®
Ultra-inert
UHPLC & HPLC Columns


Application #**AN1610**

Conditions

Column:	ACE Excel 3 C18-PFP
Dimensions:	150 x 4.6 mm
Part Number:	EXL-1110-1546U
Mobile Phase:	A: H_2O
	B: MeCN
	C: 10% formic acid

Time (mins)	%A	%B	%C
0	84	15	1
12	74	25	1
14	59	40	1
16	84	15	1
18	84	15	1

Flow Rate:	1 mL/min
Detection:	UV, 268 nm



1. Sulfadiazine
2. Sulfapyridine
3. Sulfamerazine
4. Sulfamoxole
5. Sulfamethazine
6. Sulfamonomethoxine
7. Sulfachloropyridazine
8. Sulfamethoxazole
9. Sulfadimethoxine

Reproduced with permission of National Research Institute of Agriculture and Food Technology INIA, Spain

© 2015 ACT Ltd. All rights reserved. All products are available worldwide.