
USP Monograph - Hydroquinone

ACE®
Ultra-inert
UHPLC & HPLC Columns

Application #AN1780

Conditions

Column:	ACE Excel 2 C18
Dimensions:	50 x 3.0 mm
Part Number:	EXL-101-0503U
Mobile Phase:	H_2O/MeOH (45:55 v/v)
Flow Rate:	0.45 mL/min
Injection:	0.9 µL
Temperature:	Ambient (22 °C)
Detection:	UV, 280 nm



1. Hydroquinone

Translated method according to USP <621> guidance. Always check the latest and official method information from the relevant pharmacopoeia prior to analysis.

© 2015 ACT Ltd. All rights reserved. All products are available worldwide.