
Determination of Calycosin-7-O-β-D-glucoside and Spinosin in Shenqi Wuweizi Capsules

ACE® Ultra-inert UHPLC & HPLC Columns


Application #**AN6620**

Conditions


Columns:	ACE Excel 5 SuperC18
Dimensions:	250 x 4.6 mm
Part Numbers:	EXL-1211-2546U
Mobile Phase:	MeCN/H_2O (15:85 v/v)
Flow Rate:	1.0 mL/min
Injection:	10 µL
Temperature:	Ambient
Detection:	UV, 219 nm
Method:	ChP.2015, Vol. 1



1. Calycosin-7-O-β-D-glucoside



2. Spinosin



© 2018 ACT Ltd. All rights reserved. All products are available worldwide.